



Heraeus Med GmbH 56078119

**FITTING DETAILS**

Weight 74kg  
 Voltage  
 Frequency  
 Rating  
 PF  
 Magnetic Rating

Watertight To None  
 Cable Entry Gland  
 Method of Mounting See Figure  
 Maximum Cable Size  
 EMC  
 Lamps

**COMPONENT DETAILS**

COMPONENT	QTY	NSN	CLASS GROUP	REMARKS
TBA				

**GENERAL NOTES**

1. To be used with Medical Light Head 2003i 6515-99-524-6201.
2. Support arm can be lowered.

**HAZARDOUS ITEM** – None

**Operating theatre light rail –  
 6530-12-346-3942 - CL/GP 0854**